

**New stroke assessment panel.** The Stroke assessment Power Tool is no longer active.

## Stroke

Stroke Assessment  
Done:

No

Facial Droop :

Normal

Right

Left

Unable to  
Assess

Speech :

Normal

Abnormal

Unable to  
Assess

Arm Drift :

Normal

Right Drifts  
Down

Left Drifts  
Down

Unable to  
Assess

Visual Fields:

Normal

Abnormal

Unable To  
Comply

Cerebellar Assessment:

Normal

Abnormal

Unable to  
complete

Date/Time Last Known  
Well:




## Driver's License Reader for Elite

Using your camera requires an active incident form with the Driver's License 2D Camera Scanner button on the form template.

1. On the iPad with the iOS Elite Utilities app installed, open Elite Field.


**NOTE:** Internet is not required to use the iPad camera for scanning driver's licenses.

2. Create or open an existing incident.
3. **Locate and click** the *Driver's License 2D Camera Scanner* button.



NOTE: The button will be inactive (greyed out) if iOS Elite Utilities app is not installed

4. Select *Open*.  
The iPad switches to the camera app.



5. Fill the iPad's screen, from left to right, with the driver's license's barcode. As soon as the app reads the barcode, you are redirected to the incident form.



Click *OK* to add the imported information to the incident.

**Driver's License Scan Results**

First Name:	PAUL
Middle Name:	
Last Name:	CILLAM
Date Of Birth:	09/09/1979
Gender:	Male
Address:	108 GEM COURT
City:	City of Fort Myers
State:	FL
Postal Code:	33911
Country:	USA
Weight (lbs):	
Weight (kg):	
Driver's License Number:	F40000
State Issuing Driver's License:	FL

OK Cancel